

PM-673 ROBOCLEAN TAB

COAGULANT

PM-673 ROBOCLEAN TAB: Tablets that improve the cleaning capacity of robots. They increase the size of the finest particles allowing the robot filter media (bag or cartridge) to retain them, minimizing the cloud of particles. Chemically neuter.



METHOD OF USE AND DOSAGE

- Add a tablet (without its protective film) inside the filter of the robot. If it has two compartments break it and place a part in each of them.
- Immerse and turn on the robot inside the water.
- Recommendation: use whenever there are fine particles in the bottom.
- Important: after each cleaning cycle clean the filters deeply.

PHYSICAL PROPERTIES

- Appearance: 16 g orange tablets
- Odour: Odourless
- Water Solubility: miscible in all proportions.

PRECAUTIONARY STATEMENTS

- Keep out of reach of children.
- Do not get in eyes, on skin or on clothing.
- In case of accident or if you feel unwell seek medical advice immediately.

PACKAGING

Box of 6 units.